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Bibliografia

- 1 Mullens W, Grimm RA, Verga T, et al. Insights from a cardiac resynchronization optimization clinic as part of a heart failure disease management program. *J Am Coll Cardiol*. 3 de março de 2009;53(9):765-773.
- 2 Mittal S, Plummer C, Klepfer R, Liu S, Frank C. A Novel Algorithm Increases Effective Left Ventricular Pacing During Atrial Fibrillation in Patients Undergoing Cardiac Resynchronization Therapy: Primary Results of the Randomized CRTee Trial. Apresentado no Congresso da ESC de 2016.
- 3 Hayes DL, Boehmer JP, Day JD, et al. Cardiac resynchronization therapy and the relationship of percent biventricular pacing to symptoms and survival. *Heart Rhythm*. Setembro de 2011;8(9):1469-1475.
- 4 European Heart Rhythm Association, European Society of Cardiology, Heart Rhythm Society, et al. 2012 EHRA/HRS expert consensus statement on cardiac resynchronization therapy in heart failure: implant and follow-up recommendations and management. *Heart Rhythm*. Setembro de 2012;9(9):1524-1576.
- 5 Brignole M, Auricchio A, Baron-Esquivias G, et al. 2013 ESC Guidelines on cardiac pacing and cardiac resynchronization therapy: the Task Force on cardiac pacing and resynchronization therapy of the European Society of Cardiology (ESC). Developed in collaboration with the European Heart Rhythm Association (EHRA). *Eur Heart J*. Agosto de 2013;34(29):2281-2329.
- 6 Starling RC, Krum H, Bril S, et al. Impact of a Novel Adaptive Optimization Algorithm on 30-Day Readmissions: Evidence From the Adaptive CRT Trial. *JACC Heart Fail*. Julho de 2015;3(7):565-572.
- 7 Martin D, et al. Clinical Outcomes with Adaptive Cardiac Resynchronization Therapy: Long-term Outcomes of the Adaptive CRT Trial. Late-breaking Clinical Trial. Apresentado na HFSA 2013.
- 8 Varma N, Stadler RW, Ghosh S, Kloppe A. Influence of automatic frequent pace-timing adjustments on effective left ventricular pacing during cardiac resynchronization therapy. *Europace*. Publicado em 26 de maio de 2016.
- 9 Biffi, et al. Utilizing Short Spacing between Quadripolar LV Lead Electrodes to Avoid PNS. *Cardiostim 2014*, Sessão de Apresentação 56P. Session 56P.
- 10 Biffi M, Foerster L, Eastman W, et al. Effect of bipolar electrode spacing on phrenic nerve stimulation and left ventricular pacing thresholds: an acute canine study. *Circ Arrhythm Electrophysiol*. 1º de agosto de 2012;5(4):815-820.
- 11 Biffi M, Zanon F, Bertaglia E, et al. Short-spaced dipole for managing phrenic nerve stimulation in patients with CRT: the "phrenic nerve mapping and stimulation EP" catheter study. *Heart Rhythm*. Janeiro de 2013;10(1):39-45.
- 12 Lunati MG, Gasparini M, Landolina M, et al. Long-Term Effect of Steroid Elution on the Electrical Performance of Coronary Sinus Leads for Cardiac Resynchronization Therapy. Apresentado na HRS 2012 (AB10-05).
- 13 Demmer W. VectorExpress Performance Results. Dados da Medtronic em arquivo, janeiro de 2013.
- 14 Yarnykh VL, Terashima M, Hayes CE, et al. Multicontrast black-blood MRI of carotid arteries: comparison between 1.5 and 3 tesla magnetic field strengths. *J Magn Reson Imaging*. Maio de 2006;23(5):691-698.
- 15 Saleh RS, Lohan DG, Nael K, et al. Cardiovascular MRI at 3T. *Appl Radiol*. 12 de novembro de 2007;36(11):10-26.
- 16 Taïeb S, Devise V, Pouliquen G, et al. Clinical utility and economic viability of a 3T MRI in an anti-cancer centre: The experience of the centre Oscar Lambret. *Diagn Interv Imaging*. Julho de 2012;93(7-8):561-568.

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